COVID-19 (SARS-CoV-2) Antigen Test Kit Instructions for Use

- Self-test used on-site or at home
- Please read the instructions for use before use.

Product name

COVID-19 (SARS-CoV-2) Antigen Test Kit

Specification

1 Test/Kit

Intended use

This kit is intended for in vitro qualitative detection of the nucleocapsid protein antigen from SARS-CoV-2 in nasal swabs from individuals suspected of COVID-19. The kit is suitable for layperson use, children under age of 14 should be assisted by an adult. Please always seek for assistance and supervision with the test when needed and follow your local guidelines for specimen collection by children.

COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 2 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

Kit components

1 Antigen test cassette
1 Antigen extract R1
1 Sample swab
1 Instructions for use
1 Biohazard waste bag

Test principle

This kit employs immunochromatography for detection. The specimen will move forward along the test cassette under capillary action. If the SARS-CoV-2 viral antigen is present, it will be bounded to the colloidal gold-labeled SARS-CoV-2 specific antibodies. The immune complex will be captured by coronavirus monoclonal antibody fixed in the T (Test) line. If the specimen is SARS-CoV-2 positive, both the T (Test) line and C (Control) line will become visible. If the specimen is SARS-CoV-2 negative, the C (Control) line will become visible but the T (Test) line will be invisible. The C line must be visible if the test has been performed correctly.

Storage & stability

1. Store at 2°C-30°C. DO NOT FREEZE.
2. After the aluminum foil bag is unsealed, the test cassette should be used as soon as possible.

Preparatory steps

1. Disinfect the surface where you will open the test kit. Remove and lay out contents of the test kit on a clean, flat surface.

- Please blow your nose and clear the nostril before taking the test.
- Wash hands with soap and water. If soap and water are not available, use hand sanitizer. Dry your hand completely before taking the test.
- You need to prepare a timer or any device with timer function.

Test methods

- Press along the dotted line of the tube stand, and place the extraction tube on the tube stand. Remove the sealing film of the tube.

2. Collect the sample:

A. Remove the swab from the package. Do not touch the soft end (swab head) with your hands or anything else.

B. Insert the entire soft end of the swab into the nostril (about 1.5-2.0 cm).

C. Slowly rotate the swab, gently press against the inside of the nostril at least 5 times for a total of 15 seconds. Get as much nasal secretion as possible on the soft end of the swab.

D. Gently remove the swab.

E. Using the same swab, repeat steps B-D in the other nostril with the same end of the swab.

- Put the swab specimen into the extraction tube, rotate the swab in the liquid for about 10 seconds, and press the swab head against the tube wall to release the specimen in the swab.


- 1 -

Negative

A negative test result does not rule out possibility of infection, an infection may be present even if your test result is negative. You should continue to comply with all applicable rules regarding contact with others and protective measures. In case of suspicion, please repeat the test after 1 - 2 days, as the novel coronavirus may not be detected accurately in all phases of an infection.

• **Invalid result:** if the C line is not observed, it will be invalid regardless of whether there is T line (as shown in the figure below), and the test should be repeated.

Invalid

An invalid test result is possibly caused by faulty test execution, please repeat the test; if test result is still invalid, contact a doctor or COVID-19 test center for their professional opinions and immediately contact the manufacturer or local supplier.

Benefits and Limitations

Benefits:

- This test kit is user-friendly, does not require professional testing equipment and can be used by non-professionals on-site or at home. The results can be observed with the naked eye in only 15 minutes.
- The test result from this test can help your healthcare provider make informed recommendations for your treatment / care and help limit the spread of COVID-19 to your family and others around you.

Limitation:

- This test kit is used for in vitro diagnosis only.
- This test kit is only used to detect human anterior nasal swab extracts. The results of other specimens may be inaccurate.
- This test kit is only used for qualitative detection and cannot indicate the level of novel coronavirus antigen in the specimen.
- This test kit is only a clinical auxiliary diagnostic tool. If the result is positive, it is recommended to use other methods for further examination in time and the doctor's diagnosis shall prevail.

Clinical performance


1. **LoD:** The limit of detection (LOD, direct swab method) of COVID-19 (SARS-CoV-2) Antigen Test Kit is confirmed as 5x10³ TCID₅₀/ml. The limit of detection (LOD, direct dilution method) of COVID-19 (SARS-CoV-2) Antigen Test Kit is confirmed as 1x10³ TCID₅₀/ml.

2. Analytical specificity

1) Cross-reactivity:

No false positive test results for either COVID-19 (SARS-CoV-2) Antigen Test Kit were observed on specimens from the following disease states or specific conditions: Parainfluenza virus-1, Parainfluenza virus-2, Parainfluenza virus-3, Parainfluenza virus-4, Influenza A H1N1 virus (2009), Influenza B virus (Victoria strain), Influenza B virus (Yamagata strain), Enterovirus CA16, Respiratory syncytial virus, Mycoplasma pneumoniae, 50% nasal cleaning solution specimen of healthy population, Measles virus, Mumps virus, Coronavirus 229E, Coronavirus NL63, Coronavirus OC43, MERS-CoV, RSV-A, RSV-B.

- Squeeze the swab over the head while taking the swab out of the extraction tube to remove as much liquid as possible from the swab. Put the swab after use into the biohazard waste bag.
- Install the dropper head on the extraction tube, unseal the package and take out the antigen test cassette.
- Add 2-3 drops into the specimen well of the test cassette, and start the timer. The drops should be liquid not foams/bubbles. If the drop is foam/bubbles, add another drop.

- Read the result in 15 minutes. Strong positive results can be reported within 15 minutes, however, negative results must be reported after 15 minutes, and the results after 25 minutes are no longer valid.

- The antigen test cassette, antigen extract R1, antigen extraction tube (with dropper head) and sample swab after use should be placed in the biohazard waste bag and disposed of with household waste.

- Wash hands or re-apply hand sanitizer after taking the test.

Interpretation of test results and instructions for actions

• **Positive result:** if both the C line and the T line appear (as shown in the figure below), novel coronavirus antigen has been detected and the result is positive. Look very closely! The T line can be very faint. Any pink/purple line visible here indicates a positive result.

Here are some examples of the colors of T line.

Positive

A positive test result suggests there is currently a suspicion of COVID-19 infection, if your test result is positive, you should immediately contact a doctor/family doctor or the local health authority and ask for their professional opinions, you should comply with local self-isolation guidelines and have a COVID-19 nucleic acid PCR confirmatory test carried out to confirm the infection.

• **Negative result:** if there is only a C line (as shown in the figure below), the T line is colorless, indicating that novel coronavirus antigen has not been detected and the result is negative.

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Osmunda Medical Technology Service GmbH Treuskowallee 108, 10318 Berlin, Germany Tel: 0049-30-81865123
H&W Biotech Oy Address: Teknikantie 12, 02150, Espoo, Finland Tel: +358458833889, +3584057